


PD TEMP

High Flow Curved Stylet Catheter for Acute Peritoneal Dialysis



Specifications:

- ***Quick, easy and safe***
- ***Bent Abdomen segment helps to advance catheter in the deep pelvic area***
- ***Made of Polyamide for outstanding tissue compatibility***
- ***Precision hole punching without burrs to prevent catheter occlusion***
- ***Catheter introduction depth control***
- ***Elbow connector prevents catheter kinking***
- ***Catheter marking to show direction of distal end to facilitate insertion***
- ***Sampling port for injections during treatment***

