


PD Y Administration Set

PDS-1000LLS



Specifications:

- ***Y-type administration set with drainage tubing***
- ***For use with two dialysis solution bags equipped with luer lock input ports***
- ***Drip chambers for viewing fill and drain flow rates***
- ***Rotating Luer lock patient line connector***
- ***Drain spike for attaching to drain bag with spike port, such as Medionics' DB-301***



Medionics International Inc., 114 Anderson Avenue, Markham, Ont L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com



Made in Canada