


PEDIATRIC IPD SET

PDS-1100 PLS



- **Y-type pediatric IPD set with 150 ml vented Burette**
- **For use with dialysis solution bag equipped with luer lock input port**
- **Drip chamber for viewing drain flow rates**
- **Rotating Male Luer lock patient line connector**
- **Roller clamps with super fine flow control**
- **Drain spike for attaching to drain bag with spike port, such as Medionics' DB-301**



Medionics International Inc., 114 Anderson Avenue, Markham, Ont L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com



Made in Canada