


# PEDIATRIC IPD SET

## PDS-1200BCCH



- *Y-type pediatric IPD set with 2 – 150 ml vented Burettes*
- *For use with two dialysis solution bags equipped with luer lock input ports*
- *Drip chambers for viewing fill and drain flow rates*
- *Rotating Male Luer lock patient line connector*
- *Roller clamps with super fine flow control*
- *Drain spike for attaching to drain bag with spike port, such as Medionics' DB-301*



Medionics International Inc., 114 Anderson Avenue, Markham, Ont L6E 1A5 Canada  
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com [www.medionics.com](http://www.medionics.com)



Made in Canada