


PD Drain Bag

DB-301

Specifications:

- 3.5 Litre Drain Bag
- Input port for secure connection of the drain line spike of the "Y" administration set
- Output port with clamp to empty the drain bag
- Calibration marking for measurement of the drain volume
- Hook and string for hanging the drain bag



Medionics International Inc., 114 Anderson Avenue, Markham, Ont L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com



Made in Canada