

EASY CARE

Making APD Affordable



Very Versatile

- For APD, IPD, and TIDAL
- For adults and pediatrics
- Day or Night Exchanges

Multi-functional


- *The ideal way to Initiate early incremental dialysis !*
EASY CARE is so simple and inexpensive to use that it eliminates the technical and financial considerations in deciding whether to initiate early incremental dialysis.
- *Self-start mode to improve CAPD weekly clearances !*
EASY CARE can also be run in self-start mode to perform, very inexpensively, 1 or 2 additional night exchanges at the preset time as the residual renal function declines.
- *Prescription dialysis to suit patient's lifestyle !*
EASY CARE can be conveniently switched between day and night exchanges or any combination of them to adjust to patient's needs and lifestyle.

Very Safe & Quiet

- Gravity flow control
- No noisy pump
- No forced & pressurized painful drainage

Very Simple & Portable

- Drop in cartridge – very simple and quick set-up
- Extended Battery backup and treatment data storage
- Truly portable, weighs under 8 kg with mobile stand, ideal for travelling
- Completely mobile around the house



MEDIONICS INTERNATIONAL INC.
114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada
Tel: (905) 472-6544 info@medionics.com www.medionics.com



Made in Canada